

<b>pts panels<sup>®</sup></b> Multi-Chemistry Controls    Range Card		Level 1		Level 2	
<b>Lot #:</b>	<b>MC47</b>	<b>mg/dL</b>	<b>mmol/L</b>	<b>mg/dL</b>	<b>mmol/L</b>
<b>Expiration Date:</b>	<b>2024-07-31</b>				
<b>Single-Analyte Test Strips</b>					
Cholesterol		120-190	3.11-4.92	160-270	4.14-6.99
HDL Cholesterol* <i>(CardioChek<sup>®</sup> Plus analyzer version 1.09 or above)</i>		≤25-50	≤0.65-1.30	30-70	0.78-1.81
Triglycerides		110-180	1.24-2.03	160-270	1.81-3.05
Glucose/eGLU		20-110	1.11-6.11	190-330	10.55-18.32
<b>Multi-Analyte Test Strips</b>					
Cholesterol		100-240	2.59-6.22	160-330	4.14-8.55
HDL Cholesterol* <i>(CardioChek<sup>®</sup> Plus analyzer version 1.09 or above)</i>		≤15-50	≤0.39-1.30	25-80	0.65-2.07
Triglycerides		80-160	0.90-1.81	150-270	1.70-3.05
Glucose		30-115	1.67-6.38	110-340	6.11-18.87
* HDL Cholesterol requires separate HDL Cholesterol Controls (REF 722) for CardioChek Plus (version 1.08 or below), CardioChek PA, and CardioChek analyzers.					
PTS Diagnostics Customer Service: 877-870-5610   Fax: 317-870-5608   Outside of US: +1-317-870-5610 Email: customerservice@ptsdiagnostics.com					
PS-002519 Rev. 15					