

pts panels[®]

Multi-Chemistry Controls Range Card

Level 1

Level 2

Lot #: MC48
Expiration Date: 2025-02-28

mg/dL

mmol/L

mg/dL

mmol/L

Single-Analyte Test Strips

Cholesterol	120-200	3.11-5.18	150-300	3.89-7.77
HDL Cholesterol* <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>	≤25-50	≤0.65-1.30	30-70	0.78-1.81
Triglycerides	100-180	1.13-2.03	160-260	1.81-2.94
Glucose/eGLU	40-110	2.22-6.11	150-340	8.33-18.87

Multi-Analyte Test Strips

Cholesterol	105-230	2.72-5.96	140-330	3.63-8.55
HDL Cholesterol* <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>	≤15-50	≤0.39-1.30	25-80	0.65-2.07
Triglycerides	80-160	0.90-1.81	150-260	1.70-2.94
Glucose	30-115	1.67-6.38	140-340	7.77-18.87

* HDL Cholesterol requires separate HDL Cholesterol Controls (REF 722) for CardioChek Plus (version 1.08 or below), CardioChek PA, and CardioChek analyzers.

PTS Diagnostics Customer Service: 877-870-5610 | Fax: 317-870-5608 | Outside of US: +1-317-870-5610
Email: customerservice@ptsdiagnostics.com

PS-002519 Rev. 15