

pts panels[®] Multi-Chemistry Controls Range Card		Level 1		Level 2	
		mg/dL	mmol/L	mg/dL	mmol/L
Lot #:	MC49				
Expiration Date:	2025-04-30				
Single-Analyte Test Strips					
Cholesterol		120-200	3.11-5.18	140-290	3.63-7.51
HDL Cholesterol* <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>		≤25-50	≤0.65-1.30	30-65	0.78-1.68
Triglycerides		100-180	1.13-2.03	160-260	1.81-2.94
Glucose/eGLU		40-110	2.22-6.11	150-380	8.33-21.09
Multi-Analyte Test Strips					
Cholesterol		90-220	2.33-5.70	140-320	3.63-8.29
HDL Cholesterol* <i>(CardioChek[®] Plus analyzer version 1.09 or above)</i>		≤15-60	≤0.39-1.55	25-75	0.65-1.94
Triglycerides		80-160	0.90-1.81	150-260	1.70-2.94
Glucose		40-115	2.22-6.38	150-330	8.33-18.32
<p>* HDL Cholesterol requires separate HDL Cholesterol Controls (REF 722) for CardioChek Plus (version 1.08 or below), CardioChek PA, and CardioChek analyzers.</p> <p>PTS Diagnostics Customer Service: 877-870-5610 Fax: 317-870-5608 Outside of US: +1-317-870-5610 Email: customerservice@ptsdiagnostics.com</p>					