

pts panels[®]

Multi-Chemistry Controls

Range Card

Level 1

Level 2

Lot #:

MC51

Expiration Date:

2025-08-31

mg/dL

mmol/L

mg/dL

mmol/L

Single-Analyte Test Strips

Cholesterol

110-220

2.85-5.70

150-280

3.89-7.25

HDL Cholesterol*

(CardioChek[®] Plus analyzer version 1.09 or above)

≤25-50

≤0.65-1.30

30-70

0.78-1.81

Triglycerides

100-180

1.13-2.03

170-270

1.92-3.05

Glucose/eGLU

30-110

1.67-6.11

160-380

8.88-21.09

Multi-Analyte Test Strips

Cholesterol

100-220

2.59-5.70

130-300

3.37-7.77

HDL Cholesterol*

(CardioChek[®] Plus analyzer version 1.09 or above)

≤15-50

≤0.39-1.30

25-75

0.65-1.94

Triglycerides

70-170

0.79-1.92

150-270

1.70-3.05

Glucose

30-115

1.67-6.38

150-340

8.33-18.87

* HDL Cholesterol requires separate HDL Cholesterol Controls (REF 722) for CardioChek Plus (version 1.08 or below), CardioChek PA, and CardioChek analyzers.

PTS Diagnostics Customer Service: 877-870-5610 | Fax: 317-870-5608 | Outside of US: +1-317-870-5610

Email: customerservice@ptsdiagnostics.com

PS-002519 Rev. 15